infochroma ag

Sumpfstr. 3, CH-6300 Zug; Fax: 041 748 50 65, Tel: 041 748 50 60
e-mail: info@infochroma.ch, www.infochroma.ch

Sample filtration for UPLC

The smaller the particle size of the HPLC phase the more important is the purity of the sample. In the „good old times“ of 5 μm phases, filtration wasn't always necessary. But today, with modern UPLC and its particle sizes of 1.7 μm , filtering the sample with a filter of 0.2 μm pore size is essential in order to give the column a reasonable life span.

A PTFE membrane is recommended for organic solvents, Regenerated Cellulose for aqueous protein samples. The chosen pore size should be 0.2 μm so that viruses, bacteria and other biological fragments are effectively retained. The necessary filter diameter is defined by the volume of the sample and the load of suspended particles. Various filter diameters are available so that the optimum filter can be chosen.

For applying the sample on the filter, PP-syringes with luer-slip with volumes from 1 to 50 ml are available depending on filtration volume. We are happy to send you samples of filter and matching PP-syringe for your filtrate. Write to us for samples.

UPLC/Nano HPLC is also a topic at the up-coming infoExpo in Lausanne on 23. September 2008 and Basel on 24. September 2008. A good chance to learn more about the latest developments on this subject in the seminar and on the practical aspects in the afternoon workshops.

We look forward to seeing you either in Lausanne or in Basel and to discussing your applications and possible solutions with you.

Filter with PTFE Membrane for organic samples Pore Size: 0.2 μm

≤ 2 ml Sample Volume	≤ 10 ml Sample Volume	≤ 20 ml Sample Volume	≤ 100 ml Sample Volume	≤ 200 ml Sample Volume
				
4 mm \varnothing	13 mm \varnothing	17 mm \varnothing	25 mm \varnothing	30 mm \varnothing
8804C-P-2	8813C-P-2	8817-P-2	8825C-P-2	8830-P-2
CHF 107.00/100 pcs.	CHF 148.00/100 pcs.	CHF 330.00/200 pcs.	CHF 162.00/100 pcs.	CHF 230.00/100 pcs.

Filter with Regenerated Cellulose Membrane for aqueous, protein samples Pore Size: 0.2 μm

≤ 2 ml Sample Volume	≤ 10 ml Sample Volume	≤ 20 ml Sample Volume	≤ 100 ml Sample Volume	≤ 200 ml Sample Volume
				
4 mm \varnothing	13 mm \varnothing	17 mm \varnothing	25 mm \varnothing	30 mm \varnothing
8804C-RC-2	8813C-RC-2	8817-RC-2	8825C-RC-2	8830-RC-2
CHF 107.00/100 pcs.	CHF 155.00/100 pcs.	CHF 330.00/200 pcs.	CHF 186.00/100 pcs.	CHF 220.00/100 pcs.

Polypropylene Syringes with Luer Slip

1 ml with Luer Slip 	8899-PP1 = CHF 57.00/100 pcs. <i>Special price* = CHF 44.00/100 pcs.</i>	20 ml with Luer Slip 	8899-PP20 CHF 75.00/100 pcs. <i>Special price* CHF 57.00/100 pcs.</i>
3 ml with Luer Slip 	8899-PP3 = CHF 38.00/100 pcs. <i>Special price* = CHF 29.00/100 pcs.</i>	50 ml with Luer Slip 	8899-PP50 CHF 170.00/100 pcs. <i>Special price* CHF 116.00/100 pcs.</i>
5 ml with Luer Slip 	8899-PP5 = CHF 35.00/100 pcs. <i>Special price* = CHF 26.50/100 pcs.</i>		
10 ml with Luer Slip 	8899-PP10 = CHF 56.00/100 pcs. <i>Special price* = CHF 43.00/100 pcs.</i>		

**Special price valid if filter and syringes are purchased in combination
valid only if order code **FS1024** is quoted
offer valid till 24. October 2008*

Filter: Price for 1 Pack; Quantity breaks from 3 packs
Syringes: Price for 500 pcs.; Quantity breaks from 1000 pcs.
all prices excl. VAT., excl. postage and packaging

infoExpo 2008

23. Septembre 2008

Lausanne

Hôtel Mövenpick

09:00 - 13:00: Séminaire

Detection of impurities with HPTLC

Dr. Anita Ankli, CAMAG

High-speed countercurrent chromatography (HSCCC) - an all-liquid method for the separation of natural products

Dr. Andrew Marston, Université de Genève, Laboratory Pharmacognosy and Phytochemistry

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Roland Locher, Dionex (Switzerland) AG

NQAD - Nano Quantity Analyte Detector: A new aerosol based HPLC-Detector with ultrahigh sensitivity, excellent dynamic range and good linearity

Dr. Jürgen Maier-Rosenkranz, Grace Davison Discovery Sciences Alltech Grom GmbH

MEPS, a New Possibility for Automatic Sample Preparation

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A Triple Quadrupole for GC and HPLC-MS - Technology and Application

Dr. Karl Abele, Varian AG

New techniques for accelerating throughput and efficiency in HPLC

Luca Piatti, Perkin Elmer Italy

14:15 - 17:00: Workshops

Different strategies for successful fast HPLC

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Jürg Leuenberger, CAMAG

Use and application of special phases for better separation in flash chromatography

Dr. Jörg Heerklotz / IG Instrumenten-Gesellschaft AG

Innovation in shape Selectivity HPLC for aromatic compounds and specific application for PAHs

Joyce Link, Milian SA

New techniques for accelerating throughput and efficiency in HPLC and Headspace GC

Denis Rousseil, PerkinElmer (Schweiz) AG

The practical approach to HSCCC

Ruedi Schwarz, sertec electronics

LC Techniques for Polar Compounds

Gunnar Sturm, Stagroma AG

Increasing your productivity by automating LC or GC method development

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Dr. Jürgen Maier-Rosenkranz, Grace Davison Discovery Sciences Alltech Grom GmbH

Halo. New fused core particles for very fast and very efficient HPLC separations

André Dams, Advanced Materials Technology Inc.

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Bernd Pfeffer, SGE (Deutschland) GmbH; Stand infochroma ag
Pascal Hendry, Varian AG
BGB Analytik AG
PerkinElmer (Schweiz)AG
VWR International AG