



Sumpfstr. 3, CH-6300 Zug; Fax: 041 748 50 65, Tel: 041 748 50 60  
e-mail: info@infochroma.ch, www.infochroma.ch

## GC & LC Autosampler Syringes für CTC

### GC Autosampler Syringes für CTC GC PAL ..... Page 2 - 3

- with fixed or removable needle, 0.47 mm OD, 50 mm length
- with fixed or removable needle, 0.47 mm OD, 50 mm length, multipack of 6 syringes
- with fixed or removable needle, 0.63 mm OD, 50 mm length
- with fixed or removable needle, 0.63 mm OD, 50 mm length, multipack of 6 syringes

### GC Autosampler Syringes for CTC GC PAL ..... Page 4


- with fixed or removable needle, 0.47 mm or 0.63 mm OD, various needle lengths and tips

### GC Autosampler Syringes for CTC Headspace CombiPAL ..... Page 5

- with fixed Needle 0.47 mm or 0.63 mm OD, 50 mm length

### LC Autosampler Syringes for CTC HTS, HTC & LC PAL ..... Page 6


- with fixed or removable needle, 0.72 mm OD, 50 mm length



**25%  
DISCOUNT\***  
on SGE CTC Syringes  
Contact your authorised  
distributor for more  
information  
\*until the 30th June 2009


**Order Code SGCT**

*Discount only applicable,  
if order code is quoted on order*



## GC Autosampler Syringes for CTC GC PAL

- with fixed or removable need, 0.47 mm OD, 50 mm length




Volume	Needle	length	Gauge	OD	Needle Tip	Product No.	P/N	CHF/pcs	Promotion Price until 30.06.09
0.5 $\mu$ l	removable	50 mm	26	0.47 mm	cone	8S05R-CTC-047	000 490	155.00	<b>116.25</b>
5 $\mu$ l	fixed	50 mm	26	0.47 mm	cone	8S5F-CTC-047	001 982	114.00	<b>85.50</b>
10 $\mu$ l (gastight)	fixed	50 mm	26	0.47 mm	cone	8S10F-CTC-047-GT *	002 977	80.00	<b>60.00</b>
10 $\mu$ l	fixed	50 mm	26	0.47 mm	cone	8S10F-CTC-047 *	002 980	74.00	<b>55.50</b>
10 $\mu$ l	removable	50 mm	26	0.47 mm	cone	8S10R-CTC-047	002 982	104.00	<b>78.00</b>
10 $\mu$ l (gastight)	removable	50 mm	26	0.47 mm	cone	8S10R-CTC-047-GT	002 985	109.00	<b>81.75</b>
25 $\mu$ l (gastight)	fixed	50 mm	26	0.47 mm	cone	8S25F-CTC-047-GT	003 700	104.00	<b>78.00</b>
25 $\mu$ l	fixed	50 mm	26	0.47 mm	cone	8S25F-CTC-047	003 980	124.00	<b>93.00</b>
100 $\mu$ l (gastight)	removable	50 mm	26	0.47 mm	cone	8S100R-CTC-047-GT	005 333	145.00	<b>108.75</b>
100 $\mu$ l (gastight)	fixed	50 mm	26	0.47 mm	cone	8S100F-CTC-047-GT	005 700	119.00	<b>89.25</b>
250 $\mu$ l (gastight)	fixed	50 mm	26	0.47 mm	cone	8S250F-CTC-047-GT	006 700	117.00	<b>87.75</b>
500 $\mu$ l (gastight)	fixed	50 mm	26	0.47 mm	cone	8S500F-CTC-047-GT	007 700	123.00	<b>92.25</b>

\* auch als 6 Pack erhältlich siehe unten


- with fixed or removable need, 0.47 mm OD, 50 mm length, multipack

Volume	Needle	length	Gauge	Needle	OD	Needle Tip	P/N	CHF/pcs	Promotion Price until 30.06.09
10 $\mu$ l (gastight)	fixed	50 mm	26	0.47 mm	cone	8SK6-10F-CTC-047-GT	002 976	430.00	<b>322.50</b>
10 $\mu$ l	fixed	50 mm	26	0.47 mm	cone	8SK6-10F-CTC-047	002 986	400.00	<b>300.00</b>



## GC Autosampler Syringes for CTC GC PAL

- with fixed or removable need, 0.63 mm OD, 50 mm length




Volume	Needle	length	Gauge	OD	Needle Tip	Produkte Nr	P/N	CHF/pcs	Promotion Price until 30.06.09
0.5 $\mu$ l	removable	50 mm	23	0.63 mm	cone	8S05R-CTC-063	000 492	155.00	<b>116.25</b>
5 $\mu$ l	fixed	50 mm	23	0.63 mm	cone	8S5F-CTC-063	001 981	114.00	<b>85.50</b>
5 $\mu$ l	removable	50 mm	23	0.63 mm	cone	8S5R-CTC-063	001 984	123.00	<b>92.25</b>
10 $\mu$ l	fixed	50 mm	23	0.63 mm	cone	8S10F-CTC-063 *	002 981	74.00	<b>55.50</b>
10 $\mu$ l	removable	50 mm	23	0.63 mm	cone	8S10R-CTC-063	002 984	104.00	<b>78.00</b>
10 $\mu$ l (gastight)	fixed	50 mm	23	0.63 mm	cone	8S10F-CTC-063-GT *	002 987	88.00	<b>66.00</b>
25 $\mu$ l (gastight)	fixed	50 mm	23	0.63 mm	cone	8S25F-CTC-063-GT	003 987	123.00	<b>92.25</b>
100 $\mu$ l (gastight)	fixed	50 mm	23	0.63 mm	cone	8S100F-CTC-063-GT	005 335	145.00	<b>108.75</b>

\* auch als 6 Pack erhältlich siehe unten

- with fixed or removable need, 0.63 mm OD, 50 mm length, multipack


Volume	Needle	length	Gauge	Needle	OD	Needle Tip	6 pack	CHF/pcs	Promotion Price until 30.06.09
10 $\mu$ l	fixed	50 mm	23	0.63 mm	cone	8SK6-10F-CTC-063	002 971	400.00	<b>300.00</b>
10 $\mu$ l (gastight)	fixed	50 mm	23	0.63 mm	cone	8SK6-10F-CTC-063-GT	002 972	430.00	<b>322.50</b>

Promotion Prices = net price, no further discounts applicable. Only valid if promotion code SGCT is quoted on order.  
Prices excl. VAT, excl. postage and packaging CHF 15.00; Offer expires 30.06.2009



## GC Autosampler Syringes for CTC GC PAL


- with fixed or removable need, 0.47 mm or 0.63 mm OD, various needle lengths and tips



Volume	Needle	length	Gauge	OD	Needle Tip	Product No.	P/N	CHF/pcs	Promotion Price until 30.06.09
10 $\mu$ l	fixed	50 mm	26	0.47 mm	bevel	8S10F-CTC-047BV	002 705	90.00	<b>67.50</b>
10 $\mu$ l	fixed	70 mm	26	0.47 mm	cone	8S10F-CTC-70/047	002 978	50.00	<b>37.50</b>
10 $\mu$ l	fixed	80 mm	26	0.47 mm	cone	8S10F-CTC-80/047	002 992	88.00	<b>66.00</b>
10 $\mu$ l	removable	80 mm	26	0.47 mm	cone	8S10R-CTC-80/047	002 995	104.00	<b>78.00</b>
25 $\mu$ l (gastight)	removable	50 mm	26	0.47 mm	side hole / dome tip	8S25R-CTC-047SH-GT	003 988	166.00	<b>124.50</b>

Volume	Needle	length	Gauge	OD	Needle Tip	Product No.	P/N	CHF/pcs	Promotion Price until 30.06.09
10 $\mu$ l	fixed	80 mm	23	0.63 mm	cone	8S10F-CTC-80/063	002 989	86.00	<b>64.50</b>
100 $\mu$ l	removable	50 mm	23	0.63 mm	side hole / dome tip	8S100R-CTC-063SH-GT	005 337	169.00	<b>126.75</b>

Promotion Prices = net price, no further discounts applicable. Only valid if promotion code SGCT is quoted on order.  
Prices excl. VAT, excl. postage and packaging CHF 15.00; Offer expires 30.06.2009




## GC Autosampler Syringes for CTC Headspace CombiPAL

- with fixed needle, 0.47 mm or 0.63 mm OD, 50 mm length

Volume	Needle	length	Gauge	OD	Needle Tip	Product No.	P/N	CHF/pcs	Promotion Price until 30.06.09
1 ml	fixed	50 mm	23	0.63 mm	side hole / dome tip	8S1MF-CTC-063SH-GT	008 130	260.00	<b>195.00</b>
1 ml	fixed	50 mm	26	0.47 mm	side hole / dome tip	8S1MF-CTC-047SH-GT	008 135	290.00	<b>217.50</b>
2.5 ml	fixed	50 mm	23	0.63 mm	side hole / dome tip	8S2.5MF-CTC-063SH-GT	008 630	260.00	<b>195.00</b>
2.5 ml	fixed	50 mm	26	0.63 mm	side hole / dome tip	8S2.5MF-CTC-047SH-GT	008 635	310.00	<b>232.50</b>

*Promotion Prices = net price, no further discounts applicable. Only valid if promotion code SGCT is quoted on order.  
Prices excl. VAT, excl. postage and packaging CHF 15.00; Offer expires 30.06.2009*



## LC Autosampler Syringes for CTC HTS, HTC & LC PAL

- with fixed or removable need, 0.72 mm OD, 50 mm length

Volume	Needle	length	Gauge	OD	Needle Tip	Product No.	P/N	CHF/pcs	Promotion Price until 30.06.09
10 $\mu$ l	fixed	50 mm	22	0.72 mm	LC tip	8S10F-CTC-072LC	002 710	260.00	195.00
10 $\mu$ l (gastight)	fixed	50 mm	22	0.72 mm	LC tip	8S10F-CTC-072LC-GT	002 965	290.00	217.50
10 $\mu$ l	fixed	50 mm	22	0.72 mm	side hole / dome tip	8S10F-CTC-072SH-GT	002 983	260.00	195.00
25 $\mu$ l	fixed	50 mm	22	0.72 mm	LC tip	8S25F-CTC-072LC	003 715	310.00	232.50
25 $\mu$ l (gastight)	removable	50 mm	22	0.72 mm	LC tip	8S25R-CTC-072LC-GT	003 985	140.00	105.00
25 $\mu$ l (gastight)	removable	50 mm	22	0.72 mm	side hole / dome tip	8S25R-CTC-072SH-GT	003 989	169.00	126.75
50 $\mu$ l (gastight)	fixed	50 mm	22	0.72 mm	LC tip	8S25F-CTC-072LC-GT	004 810	174.00	130.50
100 $\mu$ l (gastight)	fixed	50 mm	22	0.72 mm	LC tip	8S100F-CTC-072LC-GT	005 715	119.00	89.25

Promotion Prices = net price, no further discounts applicable. Only valid if promotion code SGCT is quoted on order.  
Prices excl. VAT, excl. postage and packaging CHF 15.00; Offer expires 30.06.2009