

## Black-and-White or colourful?

The well known and well liked Sample Vials are used for storing samples of any kind. Thanks to our own production and a large order, we have been able to further develop this product line and concept.

So far our Sample Vial line followed American practice and standards. However, in Europe, metric measurements instead of ounces are in use and packaging is by the hundred rather than by the dozen. There were also frequent requests for colour coded caps but unfortunately we were not able to provide caps in reasonably small quantities.



Thanks to our own production, we have now invested in the necessary tools and forms and are therefore able to offer the first "sample vial family" with a diameter of 27 mm and volumes of 10, 15, 20 and 25 ml. The solid screw caps are available in 6 different colours and have a PTFE liner to keep your sample in an inert environment. Packed in high quality cardboard boxes for 100 vials, the vials are separated by a grit divider which keeps them from knocking together.

Order some samples of these storage vials or profit from our "Sample Kit" G075X-27Kit25-D and always have the proper vial at your disposal.

[http://www.infochroma.ch/eng/storagevials/storage\\_q27.html](http://www.infochroma.ch/eng/storagevials/storage_q27.html)

## Black-and-White or colourful?



Sample Kit G075X-27Kit25-D

CHF 114.00/Pack

The Kit of our first Storage Vial family:

contains 25 Stk. each Storage Vial, amber glass, ø 27 mm;  
 10 ml volume, 37 mm high  
 15 ml volume, 47 mm high  
 20 ml volume, 57 mm high  
 25 ml volume, 68 mm high

incl. closed top screw cap with PTFE liner, 6 different colours  
 - assorted

packed in our white high quality cardboard box with grit  
 divider for 100 vials




Cap Colour	Volume	Height	Price / 100 pcs.
Y	10 ml	37 mm	G075X-27/037-D CHF 114.00
B	15 ml	47 mm	G075X-27/047-D CHF 117.00
G	20 ml	57 mm	G075X-27/057-D CHF 117.00
R	25 ml	68 mm	G075X-27/068-D CHF 119.00
S	10 ml	37 mm	G075*-27/037-D CHF 119.00
W	15 ml	47 mm	G075*-27/047-D CHF 122.00
	20 ml	57 mm	G075*-27/057-D CHF 122.00
	25 ml	68 mm	G075*-27/068-D CHF 124.00

all prices are net in CHF; excl. VAT; excl. CHF 12.00 postage and packaging

### APPLICA 2006

to good customers  
 we hand out a voucher  
 for a free entry -  
 value CHF 390.00

conditions apply  
 as long as supply lasts



Dienstag  
 12. September 2006  
 Basel  
 UBS Ausbildung- und  
 Konferenzzentrum



Mercredi  
 13 septembre 2006  
 Lausanne  
 Hôtel Mövenpick

apply under [www.infomac.ch/applca06](http://www.infomac.ch/applca06)  
 more information under [www.sach.ch/applca](http://www.sach.ch/applca)

register under [www.infoexpo.ch/anmeldung](http://www.infoexpo.ch/anmeldung)  
 more information under [www.infoexpo.ch](http://www.infoexpo.ch)