

don't go green and yellow with envy....

We have extended our storage vial family by adding screw caps for 2 ml standard vials with the "Narrow Mouth" screw thread 8-425. We are now able to provide vials and caps from our own production and thus can guaranty that both go perfectly together.

At the moment, we have yellow and green closed top screw caps in stock, but other colours will follow soon. The 2 ml vials and our "Vase-Vial" with a residual volume of less than 15 μ l complete the smaller end of our storage vial family and are ideally suited for the storage of smallest sample volumes. Of course open top screw caps are also available even open top screw caps with a wider flange for autosampler where the vial hangs in the vial tray.



It is our pleasure to send you some samples. If you are looking for a special cap colour do not hesitate to contact us for the colour of your choice.

Are you daily handling small sample volumes? If so we recommend the lecture of Mr. Dr. Norbert Hell from TeLa GmbH Bremerhave on „Mikroprobenaufbereitung und Dosierung“ at the applica 2006 in Olten on Wednesday 8. November 2006. Apply for the last voucher for 8. November 2006 which gives you a free entry to applica. Please be aware that the main language of the seminar is German.

for further information please refer to

<http://www.infochroma.ch/eng/storagevials/index.html>
http://www.infochroma.ch/eng/storagevials/storage_g12.html
http://www.infochroma.ch/eng/storagevials/storage_cap.html
<http://www.infomac.ch/applica06>
<http://www.sach.ch/applica.html>



Storage vials for smallest sample volumes

- storage vials with „Narrow Mouth“ screw thread 8-425
- yellow and green closed top caps with PTFE liner; other colours will follow
- 100 % compatibility between vial and cap since both are from own production
- choose between 1.5 ml V-Vial with a residual volume smaller than 15 μ l or 2 ml Vial
- also available: open top caps or open top caps with wider flange for Autosamplers where the vial hangs in the vial tray.



www.infochroma.ch/eng/storagevials/storage_g12.html



X = cap colour according to stock

***** = replace with code for cap colour

Price / 100 pcs.

Price / 100 pcs.

Price / 100 pcs.

Price / 100 pcs.

G075**X**-12/032-H/V15 μ
EUR 69.00

G075**X**-12/032-D/V15 μ
EUR 72.00

G075**X**-12/032-H
EUR 40.00

G075**X**-12/032-D
EUR 42.00

G075***-**12/032-H/V15 μ
EUR 73.00

G075***-**12/032-D/V15 μ
EUR 76.00

G075***-**12/032-H
EUR 44.00

G075***-**12/032-D
EUR 45.00

all prices are net in EUR; quantity breaks applicable; excl. VAT; excl. EUR 10.00 postage and packaging

APPLICA 2006

to good customers
we hand out a voucher
for a free entry -
value CHF 390.00

conditions apply
as long as supply lasts

apply under www.infomac.ch/applica06
more information under www.sach.ch/applica

