

Visiting our booth at Analytica 2006 or Achema 2006 is worthwhile. In addition to our comprehensive range of autosampler vials, sample vials and HPLC filters, we introduce:

NEW: 250 μ l Vial - iV2 μ
- suitable for volatile samples -

Like the μ -Vial, the new Vial iV2 μ has the insert fused into the body of a 2 ml vial. However unlike the μ -Vial, it is the insert of the iV2 μ that is closed with the cap. That makes the iV2 μ suitable for use with volatile samples as no gas can evaporate into the space between vial and insert.



The 250 μ l insert of the iV2 μ has a very fine tip and the dead volume is less than 5 μ l. The iV2 μ can also be used on Waters Autosamplers which use syringe needles with side hole.



Visit us in
Hall 5.1
Booth A15



Analytica 2006

Visit us in Hall A2, Booth 177

Neue Messe München
25. – 28. April

NEW: FIamp
- reclosable ampule -

FIamp is filled and sealed like any conventional ampule; the sample is hermetically sealed, ideal for storage.

However after breaking the ampule, it can be reclosed with a plug. The sample does not have to be filled into another vial to be kept in a clean environment.



GTG - your partner for customized products:

In the glass trade centre of Thuringen, we develop and produce innovative quality glass products for Chromatography, Pharmacy and Industry for the German and foreign market.

We develop products according to your specifications.

International Distributors

We wish to enlarge our existing international dealer network for our own products. We name one distributor per country and are looking for distributors for those countries that are not yet represented. You find a list of our distributors on: www.gtgvials.de „Contact“

NEW: Barcode Labels for the Lab
- require exact analysis of your needs -

Barcode labels cannot be sold like a mass product. Barcode labels have to fit the requirements and routine of a laboratory. This requires an exact analysis of requirements and conditions.

Our experience with the production of barcode labels allows us to offer an exact analysis of your requirements together with a quality controlled production.



Optimize Technologies
- innovative HPLC products -

Our co-exhibitor Optimize Technologies produces innovative HPLC products, which facilitate handling in the laboratory, eg:

OPTI-LINX quick connect hardware - a high-pressure, low-dispersion quick-connect system for HPLC

OPTI-MAX
the all-in-one check valve cartridge

