



Sumpfstr. 3, CH-6300 Zug; Fax: 041 748 50 65, Tel:041 748 50 60; info@infochroma.ch

Large gob - small bottle

Ever had to take a sample and looked for the suitable bottle? But where to find one?

There is a complete range of sample bottles for all kinds of samples starting from 6 ml. Made from moulded glass, these bottles come with an inert green cap with a PTFE liner so that a sample can be stored without problems until it's analysed.

There are sample bottles with a narrow screw neck opening or "jam jar" like sample bottles with a wide opening ideal for solid samples. For the analysis of volatile components, there are open-top screw caps which allow the direct analysis of the sample.

Volumes range from 6 ml up to 4000 ml and the bottles are available in clear or amber glass. In contrast to vials produced in Gräfenroda, which are made from tubular glass, sample bottles are formed on an Individual Section (IS) Machine, which shapes the gob (a lump of molten glass) into a glass container. Therefore the walls of a sample bottles is much thicker which also makes it more robust.

Are you looking for a specific bottle? Contact us for more information or samples. Or register for the next infoExpo in Bern on Tuesday 20. March 2007 or in Zurich on Wednesday 21. March 2007 and after the interesting seminar and a delicious lunch buffet, visit our booth. We are happy to show you the various bottles and to discuss the advantages and disadvantages of the various types.



Sample Bottles & Closed-Top Screw Caps with PTFE Liner



	Boston Round Clear & Amber Glass	Wide Mouth Packer Amber Glass	Standard Wide Mouth Clear Glass	French Square Clear Glass	Graduated Clear Glass	Straight Sided Rounds Clear & Amber Glass	
6 ml		24 x 43 mm					
15 ml		31 x 50 mm		25 x 61 mm			
30 ml	31 x 79 mm*	37 x 65 mm		31 x 72 mm 34 x 61 mm	34 x 68 mm	43 x 43 mm*	
60 ml	39 x 94 mm	45 x 75 mm		38 x 86 mm 39 x 89 mm	42 x 83 mm	55 x 55 mm* 54 x 51 mm**	
120 ml	49 x 114 mm	54 x 95 mm	52 x 96 mm	44 x 144 mm	51 x 102 mm	60 x 60 mm*	
125 ml						61 x 67 mm**	
180 ml						65 x 79 mm*	
240 ml	61 x 138 mm		62 x 104 mm	56 x 137 mm	62 x 127 mm	74 x 89 mm* 70 x 92 mm**	
250 ml		65 x 120 mm					
300 ml		65 x 100 mm					
480 ml	70 x 160 mm	70 x 140 mm	72 x 133 mm	68 x 167 mm	70 x 130 mm	90 x 95 mm*	
500 ml		80 x 146 mm				108 x 89 mm**	
960 ml	90 x 190 mm	90 x 160 mm	94 x 177 mm	83 x 203 mm		90 x 160 mm*	
1000 ml	92 x 205 mm**	99 x 178 mm					
1250 ml		106 x 190 mm					
1920 ml			122 x 214 mm				
2500 ml		140 x 239 mm					
	* clear glass only ** amber glass only	www.infochroma.ch/eng/samplebottles					* clear glass only ** amber glass only



Dienstag
20. März 2007

Bern
Hotel Ambassador

Register for infoExpo 2007

Programme and Registration
www.infoexpo.ch



Mittwoch
21. März 2007

Zürich
ConventionPoint
SWX Swiss Exchange