

infochroma ag

Sumpfstr. 3, CH-6300 Zug; Fax: 041 748 50 65, Tel: 041 748 50 60
e-mail: info@infochroma.ch, www.infochroma.ch

Impure drugs are a danger to our kids, impure septa to your analysis

At this year's infoExpo you can look forward to a top-class lecture from the analytical front. Behind the title „HPLC on the dance floor“ is the mobile analytical laboratory that is present at Parties and the Streetparade and allows partygoers to have their party drugs analysed for contents and impurities. In the most harmless case it may be a pill of pure Acetylsalicylic acid, better known as Aspirin, for which 15CHF have been paid. But what if the hip drug has been adulterated with an unknown substance or one that is even more noxious than the drug itself?

We are convinced that you will get interesting answers not just to this question but also to the analysis on the road, which is what makes our profession so interesting. With your quick registration you are securing yourself a place at the infoExpo on 4th November 2009 in Basel.

Just as it is with Drugs, impurities in your Autosampler Vial Septa can be the reason for distorted results. After a year of extensive field testing with the most critical clientel, our Swiss customers, we are so happy with our ms-Pure Silicone/PTFE Septa from the infochroma production that we don't just want to advertise it, we recommend it to you. The simple principle is non-pigmented Silicone and non-pigmented PTFE for an absolutely neutral behaviour and the best possible chromatographic purity. As these septa come straight from the infochroma production with no intermediate trade, we can offer them at an absolutely competitive price. Test our ms-Pure caps and convince yourself of the excellent price-value ratio. We also took the packaging into consideration. Our ms-Pure caps come in clear PET screw-top jars. So that you can concentrate on your analysis.

We look forward to showing you our ms-Pure caps at the infoExpo in Basel on 4. November 2009 after the lecture of „HPLC on the Dance Floor“ and a tasty lunch buffet. If you can't make it to infoExpo, we are happy to send you samples or to show you the new caps in your laboratory. Contact us.

<http://www.eve-rave.net/abfahrer/drugchecking.sp?text=240&page=1>

<http://www.saferparty.ch/>

http://www.nzz.ch/nachrichten/schweiz/ein_aspirin_fuer_15_franken_an_der_street_parade_1.3302945.html

<http://www.infoexpo.ch>

http://www.infochroma.ch/eng/mspure_septum/index.html



Sumpfstr. 3, CH-6300 Zug; Fax: 041 748 50 65, Tel: 041 748 50 60
e-mail: info@infochroma.ch, www.infochroma.ch



on Wednesday 04. November in Basel

Visit us!

register now under
www.infoexpo.ch



ms-Pure Septum: Non-pigmented Silicone/PTFE Septum

With the increasing sensitivity of GC's and HPLC's there is an increased need for low bleed septum material to reduce the potential for sample error and to improve overall reliability. Pigmentation in coloured septum material can be a source of contaminant and tests show that the MS-Pure Septum is over 60 % cleaner than the standard white/red versions on the market.

ms-Pure Septum

- highest quality non-pigmented silicone elastomer, low bleed
- excellent re-sealability, resists coring
- intended for multiple injections
- low compression set for storage reliability




available as:

- Silicone/PTFE
- Silicone/PTFE pre-slit
- PTFE/Silicone/PTFE

Autosampler Vial Caps with ms-Pure Silicone/PTFE Septum

11 mm Aluminium Crimp Cap

	Septum Thickness	Septum Colour	Cap Colour	Product No.
	1 mm	clear/clear	clear green red	G003-ACC-SKFK10 G003-ACG-SKFK10 G003-ACR-SKFK10


Price at 500 pcs. = CHF 18.50/100 pcs.

11 mm Snap Cap

	Septum Thickness	Septum Colour	Cap Colour	Product No.
	1 mm	clear/clear	clear blue red	G003-SCC-SKFK10 G003-SCB-SKFK10 G003-SCR-SKFK10


Price at 500 pcs. = CHF 26.00/100 pcs.

Agilent compatible 9 mm Screw Cap

	Septum Thickness	Septum Colour	Cap Colour	Product No.
	1 mm	clear/ clear	blue red yellow black	G004-HP-CB-SKFK10 G004-HP-CR-SKFK10 G004-HP-CY-SKFK10 G004-HP-CS-SKFK10

Price at 500 pcs. = CHF 26.00/100 pcs.


8 mm Narrow Mouth Screw Cap

	Septum Thickness	Septum Colour	Cap Colour	Product No.*
	1 mm	clear/clear	yellow blue green red black white	G004-NM-CY-SKFK10 G004-NM-CB-SKFK10 G004-NM-CG-SKFK10 G004-NM-CR-SKFK10 G004-NM-CS-SKFK10 G004-NM-CW-SKFK10

Price at 500 pcs. = CHF 30.00/100 pcs.

*screw caps with flange for Japanese Autosamplers also available;
Product No = G004-NM-CJ*-SKFK10

10 mm Wide Mouth Screw Cap

	Septum Thickness	Septum Colour	Cap Colour	Product No.
	1 mm	clear/clear	black	G004-WM-CS-SKFK10

Price at 500 pcs. = CHF 36.00/100 pcs.

Price excl. VAT, excl. postage and packaging CHF 15.00
next quantity break 1000 pcs.

for further information see: www.infochroma.ch/eng/mspure_septum