



infochroma ag

Sumpfstr. 3, CH-6300 Zug; Fax: 041 748 50 65, Tel: 041 748 50 60
e-mail: info@infochroma.ch, www.infochroma.ch

Ultrafast separation of peptides with C18

PEPTIDE ES-C18 UHPLC COLUMNS



I am a sceptic, especially when it comes to American brochures and their effusive promises. I read them and, like the German comedian Mario Barth says about his girlfriend, I think: „yes, yes... don't even think about discussing it Baumli, you don't stand a chance.“ Yet, as amazing as the headlines in the new brochure about the new peptide column from HALO range sound, the boys from AMT have always kept their promises. The promise: UHPLC fast separation of peptides even with standard HPLC systems; High efficiency; Rugged and reliable; High peak capacity.

This is possible because of a bigger pore of 160Å instead of 90Å and an extra stable, steric protected, bonded C18 phase, that even at a pH below 2 and a higher temperature up to 90°C does not show any decline.

Interested? Inquire whether this column is suitable for your separation. We are happy to offer you a column for test purposes with a „refund“ warranty.

for further information please refer to
http://www.infochroma.ch/eng/hplc_columns/index.html