





**μ-Vial**

- fused insert in standard 12 x 32 mm vial
- correct vertical alignment of inserte guaranteed
- robust design
- available in clear (H) or amber (D) borosilcate glass

μ-Vial: Patent Nr. 196 12 265




**Crimp μ-Vial**




8002-CV-H/i3μ  
8002-CV-D/i3μ

**Snap/Crimp μ-Vial**




8002-SC-H/i3μ  
8002-SC-D/i3μ

**Agilent compatible Screw μ-Vial**



8004-HP-H/i3μ  
8004-HP-D/i3μ

**Wide Mouth Screw μ-Vial**



8004-WM-H/i3μ  
8004-WM-D/i3μ

Prices in CHF/Price Unit; excl. VAT		Quantity Breaks [pcs.]		
Product No.	Pack of	500 pcs	1000 pcs	5000 pcs
8002-CV-D/i3μ	500 pcs.	87.00 / 100 pcs	81.00 / 100 pcs	64.00 / 100 pcs
8002-CV-H/i3μ	500 pcs.	82.00 / 100 pcs	77.00 / 100 pcs	61.00 / 100 pcs
8002-SC-D/i3μ	500 pcs.	87.00 / 100 pcs	81.00 / 100 pcs	64.00 / 100 pcs
8002-SC-H/i3μ	500 pcs.	82.00 / 100 pcs	77.00 / 100 pcs	61.00 / 100 pcs
8004-HP-D/i3μ	500 pcs.	94.00 / 100 pcs	89.00 / 100 pcs	70.00 / 100 pcs
8004-HP-H/i3μ	500 pcs.	88.00 / 100 pcs	83.00 / 100 pcs	65.00 / 100 pcs
8004-WM-D/i3μ	500 pcs.	94.00 / 100 pcs	89.00 / 100 pcs	70.00 / 100 pcs
8004-WM-H/i3μ	500 pcs.	88.00 / 100 pcs	83.00 / 100 pcs	65.00 / 100 pcs