



infochroma ag

Sumpfstr. 3, CH-6300 Zug; Fax: 041 748 50 65, Tel: 041 748 50 60
e-mail: info@infochroma.ch, www.infochroma.ch

HPLC Syringe Filter & Accessories

Brochure 2008





infochroma ag

Sumpfstr. 3, CH-6300 Zug; Fax: 041 748 50 65, Tel: 041 748 50 60
 e-mail: info@infochroma.ch, www.infochroma.ch

Price: All Prices listed are net in CHF. They are valid at the quantity break and price unit under which they are listed. All Prices are ex works Zug excl. VAT, excl. Postage and Packaging within Switzerland and FL. CHF 13.50. Shipments abroad ex works.

Prices are subject to change without notice.

Valid: from 18.03.2008

Payment: 30 days net



samples for testing

We are happy to provide samples for test purposes. Please contact us

© Copyright:

Brochure: infochroma ag, Zug, CH
 infochroma ag, Zug, CH


Text: infochroma ag, Zug, CH

Pictures: infochroma ag, Zug, CH

Cover Page: infochroma ag, Zug, CH

Gnomes from the Gartenzwerghmanufaktur
 Philipp Griebel, Gräfenroda, D
www.zwergen-griebel.de
 in Switzerland: www.zwergshop.ch


Index	Page
Index, Terms & Conditions	2
Available Filter Diameters	3 - 4
Description of Filter Membranes	5 - 6
Order Information	7
Prices	8 - 9
HPLC Polypropylene Filter Syringes	10
Glass Syringes with Luer Lock	11



HPLC Syringe Filter and Accessories

Quality Control each batch is monitored for:


- Low UV extractables
by HPLC
- Pore Size, Burst Pressure
& Flow Rate
- Membrane Retention/
Exclusion
- Dead Volume
- tight Luer Connection










Which Filter Diameter for which Sample Volume?


Sample Volume				
≤ 2 ml	≤ 10 ml	≤ 20 ml	≤ 100 ml	≤ 200 ml
				
4 mm Filter	13 mm Filter	17 mm Filter	25 mm Filter	30 mm Filter
Sample Loss				
< 7.5 µl	< 25 µl	< 30 µl	< 50 µl	< 150 µl

HPLC Syringe Filter and Accessories




4 mm Filter		Sample Volume	Filter ø	Brand
		≤ 2 ml	4 mm	Titan or Cronus
13 & 17 mm Filter  13 mm Titan Filter no longer available 8813-N-2	 13 mm Cronus Filter ≤ 10 ml Sample Volume 8813C-N-2	≤ 10 ml	13 mm	Cronus
	 17 mm Titan Filter ≤ 20 ml Sample Volume 8817-N-2	≤ 20 ml	17 mm	Titan
25 & 30 mm Filter  25 mm Titan Filter no longer available 8825-RC-2	 25 mm Cronus Filter ≤ 100 ml Sample Volume 8825C-RC-2	≤ 100 ml	25 mm	Cronus
	 30 mm Titan Filter ≤ 200 ml Sample Volume 8830-RC-2	≤ 200 ml	30 mm	Titan

Which Pore Size?




Which Membrane?



Membranes and their properties	Tabular Summary						
	Nylon	Polyethersulfone	Polypropylene	PTFE	PVDF Poly-Vinyliden-Di-Fluoride	Regenerated Cellulose	Cellulose Acetate
Alcohol	x	x	x	x	x	x	
Aromatic Compounds				x			
Benzoide Hydrocarbons				x			
Biomolecules		x			x	x	
Solvents			x	x			
Hexane			x	x			
Organic Solvents			x	x			
Peptides					x	x	
Proteines					x	x	
weak Base	x	x	x	x	x	x	
weak Acid	x		x	x	x	x	
strong Base		x		x			
strong Acid				x			
THF	x		x	x			
Water	x	x	x	(x)	x	x	x
Aqueous Solvents	x	x	x	(x)	x	x	x
Aqueous Samples	x	x	x	(x)	x	x	x

HPLC Syringe Filter and Accessories



HPLC Syringe Filter - Description of Membranes, available Diameters and Pore Sizes

Nylon (N)

compatible with: Water, Aquous Solvents, Alcohol, Bases, THF	hydrophilic and very pure membrane, TOC analysis in trace range show an adsorption of polar substances on the membrane, especially suitable for filtration of aqueous samples, can be used up to 50°C				
incompatible with: Acids, Halogenated Solvents, Proteineous Samples with high non-specific binding affinities	4 mm	13 mm	17 mm	25 mm	30 mm
	0.2 µm	0.2 µm	0.2 µm	0.2 µm	0.2 µm
	0.45 µm	0.45 µm	0.45 µm	0.45 µm	0.45 µm
	-	-	-	-	1.2 µm
	-	-	-	-	5.0 µm

Polyethersulfone (PES)

compatible with: strong Bases, Alcohols, Proteins, Peptides	hydrophilic with very low proteine binding and high filtration rate, espacially suited for Ion Chromatography				
incompatible with: Acids, Ketone, Ester, Halogene or Benzoid Hydrocarbons	4 mm	13 mm	17 mm	25 mm	30 mm
	-	-	0.2 µm	-	0.2 µm
	-	-	0.45 µm	-	0.45 µm

Polypropylene (PP)

compatible with: Water, Aqueous Samples, Alcohols, Hexane	hydrophilic like Nylon, flowrate like Nylon, can be used up to 100°C, protein-recovery is 92 %, can be used as universal filter				
incompatible with: Bases, Acids	4 mm	13 mm	17 mm	25 mm	30 mm
	-	-	0.2 µm	-	0.2 µm
	-	-	0.45 µm	-	0.45 µm

PTFE (P)

compatible with: Solvents, Organic Solvents, Strong Acids & Bases, Aromatic Compounds	for the filtration of solvents, organic solvents, strong acids and bases, can be used up to 100°C, may only be used for aqueous samples if filter is pre-wetted, is used in TOC analysis up to 10 ml sample volume				
incompatible with: Aqueous Samples without pre-wetting of filter	4 mm	13 mm	17 mm	25 mm	30 mm
	0.2 µm	0.2 µm	0.2 µm	0.2 µm	0.2 µm
	0.45 µm	0.45 µm	0.45 µm	0.45 µm	0.45 µm
	-	-	-	-	1.0 µm

Poly-Vinyliden-Di Fluoride (PVDF)

compatible with: Alcohol, weak Acids, Proteins, Peptides and other Biomolecules	preferably used for biological samples, proteinadsorption higher than with Regenerated Cellulose, not recommended for GC-Samples or GC typical solvents				
incompatible with: Strong Acids, Bases, Ketone	4 mm	13 mm	17 mm	25 mm	30 mm
	0.2 µm	0.2 µm	0.2 µm	-	0.2 µm
	0.45 µm	0.45 µm	0.45 µm	0.45 µm	0.45 µm

Regenerated Cellulose (RC)

compatible with: Proteins, Peptides and other Biomolecules	low protein-adsorption, highly hydrophilic surface, espacially suitable for the filtration of proteins and biological samples, better chemical resistance than Cellulose-Acetate				
incompatible with: Acids, Chloroform, THF	4 mm	13 mm	17 mm	25 mm	30 mm
	0.2 µm	0.2 µm	0.2 µm	0.2 µm	0.2 µm
	0.45 µm	0.45 µm	0.45 µm	0.45 µm	0.45 µm

Cellulose Acetat (CA)

compatible with: Aqueous Samples	low protein-adsorption, hydrophilic surface				
incompatible with: Organic Solvents	4 mm	13 mm	17 mm	25 mm	30 mm
	0.2 µm	-	0.2 µm	-	0.2 µm
	0.45 µm	0.45 µm	0.45 µm	0.45 µm	0.45 µm

HPLC Syringe Filter and Accessories

4 mm, 13 mm & 17 mm HPLC Filter



Membrane	Pore Size	Colour	≤ 2 ml: 4 mm Titan Filter	≤ 2 ml: 4 mm Cronus Filter*	≤ 10 ml: 13 mm Filter*	≤ 20 ml: 17 mm Filter
Nylon (N)	0.2 μm	purple	8804-N-2	8804C-N-2	8813C-N-2	8817-N-2
	0.45 μm	green	8804-N-4	8804C-N-4	8813C-N-4	8817-N-4
Polyethersulfone (PES)	0.2 μm	dark grey	-	-	-	8817-PES-2
	0.45 μm	pale orange	-	-	-	8817-PES-4
Polypropylene (PP)	0.2 μm	dark purple	-	-	-	8817-PP-2
	0.45 μm	white	-	-	-	8817-PP-4
Teflon (PTFE)	0.2 μm	royal blue	8804-P-2	8804C-P-2	8813C-P-2	8817-P-2
	0.45 μm	yellow	8804-P-4	8804C-P-4	8813C-P-4	8817-P-4
Poly-Vinyliden-Di-Fluoride (PVDF)	0.2 μm	black	8804-PV-2	8804C-PV-2	8813C-PV-2	8817-PV-2
	0.45 μm	red	8804-PV-4	8804C-PV-4	8813C-PV-4	8817-PV-4
Regenerated Cellulose (RC)	0.2 μm	grey	8804-RC-2	-	8813C-RC-2	8817-RC-2
	0.45 μm	brown	8804-RC-4	8804-RC-4	8813C-RC-4	8817-RC-4
Cellulose Acetate (CA)	0.2 μm	pale green	8804-CA-2	-	-	8817-CA-2
	0.45 μm	orange	8804-CA-4	-	8813C-CA-4	8817-CA-4




25 mm & 30 mm HPLC Filter


* all Cronus Filters are available certified; add „Z“ to the corresponding product no., eg. 8825C-N-4/Z

Membrane	Pore Size	Colour	≤ 100 ml: 25 mm Filter*	≤ 1 μm Glass Fibre Pre-Filter	≤ 20 ml: 30 mm Filter	≤ 1 μm Glass Fibre Pre-Filter
Nylon (N)	0.2 μm	purple	8825C-N-2	-	-	8830-N-2
	0.45 μm	green	8825C-N-4	8825C-N-4/PF	-	8830-N-4
	1.2 μm	gold	-	-	8830-N-12	-
	5.0 μm	clear	-	-	8830-N-50	-
Polyethersulfone (PES)	0.2 μm	dark grey	-	-	8830-PES-2	-
	0.45 μm	pale orange	-	-	8830-PES-4	-
Polypropylene (PP)	0.2 μm	dark purple	-	-	8830-PP-2	-
	0.45 μm	white	-	-	8830-PP-4	-
Teflon (PTFE)	0.2 μm	royal blue	8825C-P-2	-	-	8830-P-2
	0.45 μm	yellow	8825C-P-4	-	-	8830-P-4
	1.0 μm	dark blue	-	-	8830-P-10	-
Poly-Vinyliden-Di-Fluoride (PVDF)	0.2 μm	black	-	-	8830-PV-2	-
	0.45 μm	red	8825C-PV-4	-	8830-PV-4	-
Regenerated Cellulose (RC)	0.2 μm	grey	8825C-RC-2	-	8830-RC-2	-
	0.45 μm	brown	8825C-RC-4	-	8830-RC-4	-
Cellulose Acetat (CA)	0.2 μm	pale green	8825C-CA-2	-	8830-CA-2	-
	0.45 μm	orange	8825C-CA-4	-	8830-CA-4	-

HPLC Syringe Filter and Accessories




Price in CHF/Price Unit (PE); excl. VAT		Quantity [Pk.]		
Product No.	Pack of	1 Pack	3 Pack	10 Pack
8804-CA-2	100 pcs.	147.00 / 1 Pk	139.00 / 1 Pk	132.00 / 1 Pk
8804-CA-4	100 pcs.	147.00 / 1 Pk	139.00 / 1 Pk	132.00 / 1 Pk
8804-N-2	100 pcs.	140.00 / 1 Pk	132.00 / 1 Pk	126.00 / 1 Pk
8804-N-4	100 pcs.	140.00 / 1 Pk	132.00 / 1 Pk	126.00 / 1 Pk
8804-P-2	100 pcs.	140.00 / 1 Pk	132.00 / 1 Pk	126.00 / 1 Pk
8804-P-4	100 pcs.	140.00 / 1 Pk	132.00 / 1 Pk	126.00 / 1 Pk
8804-PV-2	100 pcs.	147.00 / 1 Pk	139.00 / 1 Pk	132.00 / 1 Pk
8804-PV-4	100 pcs.	147.00 / 1 Pk	139.00 / 1 Pk	132.00 / 1 Pk
8804-RC-2	100 pcs.	140.00 / 1 Pk	132.00 / 1 Pk	126.00 / 1 Pk
8804-RC-4	100 pcs.	140.00 / 1 Pk	132.00 / 1 Pk	126.00 / 1 Pk
8804C-N-2	100 pcs.	107.00 / 1 Pk	103.00 / 1 Pk	100.00 / 1 Pk
8804C-N-2/Z	100 pcs.	134.00 / 1 Pk	128.00 / 1 Pk	124.00 / 1 Pk
8804C-N-4	100 pcs.	103.00 / 1 Pk	99.00 / 1 Pk	96.00 / 1 Pk
8804C-N-4/Z	100 pcs.	134.00 / 1 Pk	128.00 / 1 Pk	124.00 / 1 Pk
8804C-P-2	100 pcs.	107.00 / 1 Pk	103.00 / 1 Pk	100.00 / 1 Pk
8804C-P-2/Z	100 pcs.	134.00 / 1 Pk	128.00 / 1 Pk	124.00 / 1 Pk
8804C-P-4	100 pcs.	107.00 / 1 Pk	103.00 / 1 Pk	100.00 / 1 Pk
8804C-P-4/Z	100 pcs.	134.00 / 1 Pk	128.00 / 1 Pk	124.00 / 1 Pk
8804C-PV-2	100 pcs.	144.00 / 1 Pk	138.00 / 1 Pk	134.00 / 1 Pk
8804C-PV-2/Z	100 pcs.	179.00 / 1 Pk	172.00 / 1 Pk	167.00 / 1 Pk
8804C-PV-4	100 pcs.	144.00 / 1 Pk	138.00 / 1 Pk	134.00 / 1 Pk
8804C-PV-4/Z	100 pcs.	179.00 / 1 Pk	172.00 / 1 Pk	167.00 / 1 Pk
8804C-RC-4	100 pcs.	107.00 / 1 Pk	103.00 / 1 Pk	100.00 / 1 Pk
8804C-RC-4/Z	100 pcs.	134.00 / 1 Pk	128.00 / 1 Pk	124.00 / 1 Pk
8813C-CA-4	100 pcs.	138.00 / 1 Pk	133.00 / 1 Pk	128.00 / 1 Pk
8813C-CA-4/Z	100 pcs.	172.00 / 1 Pk	165.00 / 1 Pk	160.00 / 1 Pk
8813C-N-2	100 pcs.	131.00 / 1 Pk	126.00 / 1 Pk	122.00 / 1 Pk
8813C-N-2/Z	100 pcs.	163.00 / 1 Pk	157.00 / 1 Pk	152.00 / 1 Pk
8813C-N-4	100 pcs.	131.00 / 1 Pk	126.00 / 1 Pk	122.00 / 1 Pk
8813C-N-4/Z	100 pcs.	163.00 / 1 Pk	157.00 / 1 Pk	152.00 / 1 Pk
8813C-P-2	100 pcs.	148.00 / 1 Pk	142.00 / 1 Pk	138.00 / 1 Pk
8813C-P-2/Z	100 pcs.	185.00 / 1 Pk	177.00 / 1 Pk	172.00 / 1 Pk
8813C-P-4	100 pcs.	148.00 / 1 Pk	142.00 / 1 Pk	138.00 / 1 Pk
8813C-P-4/Z	100 pcs.	185.00 / 1 Pk	177.00 / 1 Pk	172.00 / 1 Pk
8813C-PV-2	100 pcs.	189.00 / 1 Pk	181.00 / 1 Pk	176.00 / 1 Pk
8813C-PV-2/Z	100 pcs.	240.00 / 1 Pk	226.00 / 1 Pk	219.00 / 1 Pk
8813C-PV-4	100 pcs.	172.00 / 1 Pk	165.00 / 1 Pk	160.00 / 1 Pk
8813C-PV-4/Z	100 pcs.	220.00 / 1 Pk	206.00 / 1 Pk	199.00 / 1 Pk
8813C-RC-2	100 pcs.	155.00 / 1 Pk	149.00 / 1 Pk	144.00 / 1 Pk
8813C-RC-2/Z	100 pcs.	193.00 / 1 Pk	185.00 / 1 Pk	180.00 / 1 Pk
8813C-RC-4	100 pcs.	155.00 / 1 Pk	149.00 / 1 Pk	144.00 / 1 Pk
8813C-RC-4/Z	100 pcs.	193.00 / 1 Pk	185.00 / 1 Pk	180.00 / 1 Pk
8817-CA-2	200 pcs.	300.00 / 1 Pk	281.00 / 1 Pk	267.00 / 1 Pk
8817-CA-4	200 pcs.	300.00 / 1 Pk	281.00 / 1 Pk	267.00 / 1 Pk
8817-N-2	200 pcs.	300.00 / 1 Pk	281.00 / 1 Pk	267.00 / 1 Pk








HPLC Syringe Filter and Accessories

Price in CHF/Price Unit (PE); excl. VAT.		Quantity [Pk.]		
8817-N-4	200 pcs.	300.00 / 1 Pk	281.00 / 1 Pk	267.00 / 1 Pk
8817-P-2	200 pcs.	330.00 / 1 Pk	310.00 / 1 Pk	290.00 / 1 Pk
8817-P-4	200 pcs.	330.00 / 1 Pk	310.00 / 1 Pk	290.00 / 1 Pk
8817-PES-2	200 pcs.	360.00 / 1 Pk	340.00 / 1 Pk	320.00 / 1 Pk
8817-PES-4	200 pcs.	360.00 / 1 Pk	340.00 / 1 Pk	320.00 / 1 Pk
8817-PP-2	200 pcs.	300.00 / 1 Pk	281.00 / 1 Pk	267.00 / 1 Pk
8817-PP-4	200 pcs.	300.00 / 1 Pk	281.00 / 1 Pk	267.00 / 1 Pk
8817-PV-2	200 pcs.	340.00 / 1 Pk	320.00 / 1 Pk	310.00 / 1 Pk
8817-PV-4	200 pcs.	340.00 / 1 Pk	320.00 / 1 Pk	310.00 / 1 Pk
8817-RC-2	200 pcs.	330.00 / 1 Pk	310.00 / 1 Pk	290.00 / 1 Pk
8817-RC4	200 pcs.	330.00 / 1 Pk	310.00 / 1 Pk	290.00 / 1 Pk
8825C-CA-4	100 pcs.	155.00 / 1 Pk	149.00 / 1 Pk	144.00 / 1 Pk
8825C-CA-4/Z	100 pcs.	193.00 / 1 Pk	185.00 / 1 Pk	180.00 / 1 Pk
8825C-N-2	100 pcs.	138.00 / 1 Pk	133.00 / 1 Pk	128.00 / 1 Pk
8825C-N-2/Z	100 pcs.	172.00 / 1 Pk	165.00 / 1 Pk	160.00 / 1 Pk
8825C-N-4	100 pcs.	138.00 / 1 Pk	133.00 / 1 Pk	128.00 / 1 Pk
8825C-N-4/PF	100 pcs.	193.00 / 1 Pk	185.00 / 1 Pk	180.00 / 1 Pk
8825C-N-4/PF/Z	100 pcs.	193.00 / 1 Pk	185.00 / 1 Pk	180.00 / 1 Pk
8825C-N-4/Z	100 pcs.	172.00 / 1 Pk	165.00 / 1 Pk	160.00 / 1 Pk
8825C-P-2	100 pcs.	162.00 / 1 Pk	155.00 / 1 Pk	151.00 / 1 Pk
8825C-P-2/Z	100 pcs.	202.00 / 1 Pk	194.00 / 1 Pk	188.00 / 1 Pk
8825C-P-4	100 pcs.	162.00 / 1 Pk	155.00 / 1 Pk	151.00 / 1 Pk
8825C-P-4/Z	100 pcs.	202.00 / 1 Pk	194.00 / 1 Pk	188.00 / 1 Pk
8825C-PV-4	100 pcs.	189.00 / 1 Pk	181.00 / 1 Pk	176.00 / 1 Pk
8825C-RC-2	100 pcs.	186.00 / 1 Pk	178.00 / 1 Pk	173.00 / 1 Pk
8825C-RC-2/Z	100 pcs.	240.00 / 1 Pk	222.00 / 1 Pk	215.00 / 1 Pk
8825C-RC-4	100 pcs.	185.00 / 1 Pk	178.00 / 1 Pk	173.00 / 1 Pk
8825C-RC-4/Z	100 pcs.	240.00 / 1 Pk	222.00 / 1 Pk	215.00 / 1 Pk
8830-CA-2	100 pcs.	220.00 / 1 Pk	205.00 / 1 Pk	195.00 / 1 Pk
8830-CA-4	100 pcs.	220.00 / 1 Pk	205.00 / 1 Pk	195.00 / 1 Pk
8830-N-12	100 pcs.	220.00 / 1 Pk	205.00 / 1 Pk	195.00 / 1 Pk
8830-N-2	100 pcs.	220.00 / 1 Pk	205.00 / 1 Pk	195.00 / 1 Pk
8830-N-4	100 pcs.	220.00 / 1 Pk	205.00 / 1 Pk	195.00 / 1 Pk
8830-N-50	100 pcs.	220.00 / 1 Pk	205.00 / 1 Pk	195.00 / 1 Pk
8830-P-10	100 pcs.	230.00 / 1 Pk	216.00 / 1 Pk	205.00 / 1 Pk
8830-P-2	100 pcs.	230.00 / 1 Pk	211.00 / 1 Pk	200.00 / 1 Pk
8830-P-4	100 pcs.	230.00 / 1 Pk	211.00 / 1 Pk	200.00 / 1 Pk
8830-PES-2	100 pcs.	240.00 / 1 Pk	222.00 / 1 Pk	211.00 / 1 Pk
8830-PES-4	100 pcs.	240.00 / 1 Pk	222.00 / 1 Pk	211.00 / 1 Pk
8830-PP-2	100 pcs.	220.00 / 1 Pk	205.00 / 1 Pk	195.00 / 1 Pk
8830-PP-4	100 pcs.	220.00 / 1 Pk	205.00 / 1 Pk	195.00 / 1 Pk
8830-PV-2	100 pcs.	240.00 / 1 Pk	220.00 / 1 Pk	209.00 / 1 Pk
8830-PV-4	100 pcs.	240.00 / 1 Pk	220.00 / 1 Pk	209.00 / 1 Pk
8830-RC-2	100 pcs.	220.00 / 1 Pk	205.00 / 1 Pk	195.00 / 1 Pk
8830-RC4	100 pcs.	220.00 / 1 Pk	205.00 / 1 Pk	195.00 / 1 Pk

HPLC Syringe Filter and Accessories



Polypropylene Syringes with Luer Slip or Luer Lock

	Volume	Product No.
	1 ml mit Luer Slip	8899-PP1
	1 ml mit Luer Lock	8899-PP1/LL
	3 ml mit Luer Slip	8899-PP3
	3 ml mit Luer Lock	8899-PP3/LL
	5 ml mit Luer Slip	8899-PP5
	5 ml mit Luer Lock	8899-PP5/LL
	10 ml mit Luer Slip	8899-PP10
	10 ml mit Luer Lock	8899-PP10/LL
	20 ml mit Luer Slip	8899-PP20
	20 ml mit Luer Lock	8899-PP20/LL

Glass Syringes with Luer Lock



Volume	Product No.
1 ml	8S1-MDF-LL-GT
2.5 ml	8S2.5-MDF-LL-GT
5 ml	8S5-MDR-LL-GT
10 ml	8S10-MDR-LL-GT
25 ml	8S25-MDR-LL-GT
50 ml	8S50-MDR-LL-GT
100 ml	8S100-MDR-LL-GT

Price in CHF/Price Unit (PE); excl. VAT		Quantity [Pk.]		
Product No.	Pack of	500 pcs	1000 pcs	5000 pcs
8899-PP1	100 pcs.	57.00 / 100 pcs	54.00 / 100 pcs	52.00 / 100 pcs
8899-PP10	100 pcs.	56.00 / 100 pcs	53.00 / 100 pcs	51.00 / 100 pcs
8899-PP10/LL	1000 pcs.	77.00 / 100 pcs	73.00 / 100 pcs	70.00 / 100 pcs
8899-PP20	100 pcs.	75.00 / 100 pcs	70.00 / 100 pcs	68.00 / 100 pcs
8899-PP20/LL	1000 pcs.	97.00 / 100 pcs	92.00 / 100 pcs	88.00 / 100 pcs
8899-PP3	100 pcs.	38.00 / 100 pcs	36.00 / 100 pcs	34.00 / 100 pcs
8899-PP3/LL	100 pcs.	53.00 / 100 pcs	50.00 / 100 pcs	48.00 / 100 pcs
8899-PP5	100 pcs.	35.00 / 100 pcs	33.00 / 100 pcs	31.00 / 100 pcs
8899-PP5/LL	100 pcs.	91.00 / 100 pcs	85.00 / 100 pcs	82.00 / 100 pcs
8S01-MDF-LL-GT	1 pcs.	81.00 / 1 Pk	78.00 / 1 Pk	77.00 / 1 Pk
8S025-MDF-LL-GT	1 pcs.	81.00 / 1 Pk	78.00 / 1 Pk	77.00 / 1 Pk
8S05-MDR-LL-GT	1 pcs.	96.00 / 1 Pk	93.00 / 1 Pk	91.00 / 1 Pk
8S10-MDR-LL-GT	1 pcs.	104.00 / 1 Pk	101.00 / 1 Pk	99.00 / 1 Pk
8S100-MDR-LL-GT	1 pcs.	380.00 / 1 Pk	370.00 / 1 Pk	360.00 / 1 Pk
8S25-MDR-LL-GT	1 pcs.	240.00 / 1 Pk	230.00 / 1 Pk	226.00 / 1 Pk
8S50-MR-LL-GT	1 pcs.	270.00 / 1 Pk	259.00 / 1 Pk	254.00 / 1 Pk

* Quantity Breaks 1 Pack, 3 Pack, 10 Packs