

HALO® Selectivity-Kit

Compare the selectivity of
3 HALO® phases at an unmatched price

Efficiency, retention & selectivity are the three most important parameters to influence the resolution between peaks in the chromatogram.

Modifying the selectivity using a different stationary phase while using a simple mobile phase, improves the reproducibility of a method; especially if the method is applied in different laboratories.

HALO® Selectivity Kit:

3 HALO® columns with different phases

HALO® safety glasses


Brochure with Phase description & phase comparison



90Å HALO® Small Molecule Selectivity Kit

Compare the selectivity of a HALO® Biphenyl, AQ-C18 & C18 - 2.7 µm.


Kit 1:	2.1 x 50 mm	8AM-92812-4KT	CHF	810.00/Kit
Kit 2:	2.1 x 150 mm	8AM-92812-7KT	CHF	990.00/Kit
Kit 3:	4.6 x 50 mm	8AM-92814-4KT	CHF	810.00/Kit
Kit 4:	4.6 x 150 mm	8AM-92814-7KT	CHF	990.00/Kit



1000Å HALO® Protein-Selectivity-Kit

Compare the selectivity of a HALO® Diphenyl, ES-C18 & C4 - 2.7 µm.

Kit 1:	2.1 x 50 mm	8AM-92712-4KT	CHF	960.00/Kit
Kit 2:	2.1 x 150 mm	8AM-92712-7KT	CHF	1160.00/Kit
Kit 3:	4.6 x 50 mm	8AM-92714-4KT	CHF	960.00/Kit
Kit 4:	4.6 x 150 mm	8AM-92714-7KT	CHF	1160.00/Kit



The prices in this offer are net, no further discounts applicable, excl. VAT, postage and packaging.
Offer valid until Friday 29.03.2019